




BIFURCATES



Aortic Bifurcate Prosthesis Dimensions/Sizing Guides

| Product Code | AB Size (mm) | Treatment Range (mm) | Delivery System ID (F) | Delivery System OD (F) | Ipsi Length (mm) | Contra Length (mm) |
|--------------|--------------|----------------------|------------------------|------------------------|------------------|--------------------|
| AB2298 | 22 | 17.0-19.9 | 13 | 14 | 94 | 86 |
| AB2698 | 26 | 20.0-22.9 | 13 | 14 | 94 | 86 |
| AB3098 | 30 | 23.0-26.9 | 13 | 14 | 94 | 86 |
| AB3498 | 34 | 27.0-31.0 | 15 | 16 | 94 | 86 |

Aortic Extension


| Product Code | AE Size (mm) | Treatment Range (mm) | Delivery System ID (F) | Delivery System OD (F) | Ipsi Length (mm) |
|--------------|--------------|----------------------|------------------------|------------------------|------------------|
| AE2204 | 22 | 17.0-19.9 | 13 | 14 | 42 |
| AE2604 | 26 | 20.0-22.9 | 13 | 14 | 42 |
| AE3004 | 30 | 23.0-26.9 | 13 | 14 | 42 |
| AE3404 | 34 | 27.0-31.0 | 15 | 16 | 42 |

Bifurcate and Limb Delivery Systems



The Aortic Bifurcate Delivery System working length: **54 cm**
Iliac Limb Delivery System working length: **77 cm**

LIMBS



Iliac Limb/Limb Extension Prosthesis Dimensions/Sizing Guide

| Product Code | IL Size (mm) | Treatment Range (mm) | IL Length (mm) | Delivery System OD (F) | Ipsi Length (mm) | Contra Length (mm) |
|--------------|--------------|----------------------|----------------|------------------------|------------------|--------------------|
| IL1008 | 10 | 7.0-8.9 | 82 | 12 | 128-156 | 128-147 |
| IL1010 | 10 | 7.0-8.9 | 101 | 12 | 147-175 | 147-166 |
| IL1012 | 10 | 7.0-8.9 | 120 | 12 | 166-194 | 166-185 |
| IL1014 | 10 | 7.0-8.9 | 138 | 12 | 184-212 | 184-203 |
| IL1308 | 13 | 9.0-10.9 | 82 | 12 | 128-156 | 128-147 |
| IL1310 | 13 | 9.0-10.9 | 101 | 12 | 147-175 | 147-166 |
| IL1312 | 13 | 9.0-10.9 | 120 | 12 | 166-194 | 166-185 |
| IL1314 | 13 | 9.0-10.9 | 138 | 12 | 184-212 | 184-203 |
| IL1608 | 16 | 11.0-13.9 | 82 | 12 | 128-156 | 128-147 |
| IL1610 | 16 | 11.0-13.9 | 101 | 12 | 147-175 | 147-166 |
| IL1612 | 16 | 11.0-13.9 | 120 | 12 | 166-194 | 166-185 |
| IL1614 | 16 | 11.0-13.9 | 138 | 12 | 184-212 | 184-203 |
| IL2008 | 20 | 14.0-17.9 | 82 | 12 | 128-156 | 128-147 |
| IL2010 | 20 | 14.0-17.9 | 101 | 12 | 147-175 | 147-166 |
| IL2012 | 20 | 14.0-17.9 | 120 | 12 | 166-194 | 166-185 |
| IL2014 | 20 | 14.0-17.9 | 138 | 12 | 184-212 | 184-203 |
| IL2410 | 24 | 18.0-22.0 | 101 | 13 | 147-175 | 147-166 |
| IL2412 | 24 | 18.0-22.0 | 120 | 13 | 166-194 | 166-185 |
| IL2414 | 24 | 18.0-22.0 | 138 | 13 | 184-212 | 184-203 |

For additional information and availability, please contact your Cordis representative.

For healthcare professionals only. Important information: Prior to use, refer to the "Instructions for Use" supplied with these devices for indications, contraindications, side effects, suggested procedure, warnings, and precautions. As part of the Cordis policy of continuous product development, we reserve the right to change product specifications without prior notification. For healthcare professionals only.

© Cardinal Health 2017. All right reserved. As part of the Cordis policy of continuous product development, we reserve the right to change product specifications without prior notification. The use of INCRAFT® AAA Stent-Graft System requires that physicians be specifically trained in endovascular abdominal aortic aneurysm repair techniques, including experience with high resolution fluoroscopy and radiation safety. Cordis Corporation will provide training specific to the INCRAFT® AAA Stent-Graft System. CORDIS, Cordis LOGO and INCRAFT are trademarks or registered trademarks of Cardinal Health. All rights reserved. EU530-2 05/17