

# S.M.A.R.T.™ Flex Vascular Stent System

S.M.A.R.T.™ Flex  
for your Iliac needs.

Feeling is believing.

**Cordis**  
A Cardinal Health company

S.M.A.R.T.™ Flex  
Vascular Stent System

# S.M.A.R.T.™ Flex Vascular Stent System Product Description

**Type:** Fully connected yet flexible stent design

**Material:** Nitinol, with tantalum markers on stent ends

**Maximum Guidewire:** 0.035"

**Stent Lengths:** 20 - 100 mm

**Sheath Compatibility:** 6F

**Stent Diameters:** 9-10 mm

**Stent Delivery System Working Lengths:** 80 cm (S suffix) and 120 cm (M suffix)

**Indications:** common and external iliac arteries

## Ordering Information

Product Code		Stent			Delivery System	
80cm Usable Length	120cm Usable Length	Stent Diameter (mm)	Stent Length (mm)	Recommended Lumen Size (mm)	Recommended Sheath Introducer Size	Guidewire Acceptance
SF09020SV	SF09020MV	9	20	7.5-8.5	6	.035"
SF09030SV	SF09030MV	9	30	7.5-8.5	6	.035"
SF09040SV	SF09040MV	9	40	7.5-8.5	6	.035"
SF09060SV	SF09060MV	9	60	7.5-8.5	6	.035"
SF09080SV	SF09080MV	9	80	7.5-8.5	6	.035"
SF09100SV	SF09100MV	9	100	7.5-8.5	6	.035"
SF10020SV	SF10020MV	10	20	8.5-9.5	6	.035"
SF10030SV	SF10030MV	10	30	8.5-9.5	6	.035"
SF10040SV	SF10040MV	10	40	8.5-9.5	6	.035"
SF10060SV	SF10060MV	10	60	8.5-9.5	6	.035"
SF10080SV	SF10080MV	10	80	8.5-9.5	6	.035"
	SF10100MV	10	100	8.5-9.5	6	.035"



For more information please visit [www.cardinalhealth.co.uk](http://www.cardinalhealth.co.uk)



### Important Information:

For Healthcare Professionals Only. Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information. CORDIS, Cordis LOGO and S.M.A.R.T. are trademarks of Cardinal Health and may be registered in the US and/or in other countries.

© 2020 Cardinal Health. All Rights Reserved. 100556420 03/2020