


S.M.A.R.T.™ Flex Vascular Stent System

Our inspiration came
right from the source.



Experience our Cordis
self-expanding stents.

S.M.A.R.T.™ Flex Vascular Stent System

Feeling is believing.

S.M.A.R.T.™ Flex Vascular Stent System Product Description

Type: Fully connected yet flexible stent design

Material: Nitinol, with tantalum markers on stent ends

Maximum Guidewire: 0.035"

Stent Lengths: 20-200 mm

Sheath Compatibility: 6F

Stent Diameters: 5-10 mm

Stent Delivery System Working Lengths: 80 cm (S suffix) and 120 cm (M suffix)

Ordering Information

Stent Diameter (mm)	Stent Length (mm)	Recommended Lumen Size (mm)	Recommended Sheath Introducer Size	Catalogue Number	
				Length of delivery system 80 cm	120 cm
5	30	3.5-4.5	6	SF05030SV	SF05030MV
5	40	3.5-4.5	6	SF05040SV	SF05040MV
5	60	3.5-4.5	6	SF05060SV	SF05060MV
5	80	3.5-4.5	6	SF05080SV	SF05080MV
5	100	3.5-4.5	6	SF05100SV	SF05100MV
5	120	3.5-4.5	6	SF05120SV	SF05120MV
5	150	3.5-4.5	6	SF05150SV	SF05150MV
6	30	4.5-5.5	6	SF06030SV	SF06030MV
6	40	4.5-5.5	6	SF06040SV	SF06040MV
6	60	4.5-5.5	6	SF06060SV	SF06060MV
6	80	4.5-5.5	6	SF06080SV	SF06080MV
6	100	4.5-5.5	6	SF06100SV	SF06100MV
6	120	4.5-5.5	6	SF06120SV	SF06120MV
6	150	4.5-5.5	6	SF06150SV	SF06150MV
7	30	5.5-6.5	6	SF07030SV	SF07030MV
7	40	5.5-6.5	6	SF07040SV	SF07040MV
7	60	5.5-6.5	6	SF07060SV	SF07060MV
7	80	5.5-6.5	6	SF07080SV	SF07080MV
7	100	5.5-6.5	6	SF07100SV	SF07100MV
7	120	5.5-6.5	6	SF07120SV	SF07120MV
7	150	5.5-6.5	6	SF07150SV	SF07150MV
7	200	5.5-6.5	6	SF07200SV	SF07200MV
8	30	6.5-7.5	6	SF08030SV	SF08030MV
8	40	6.5-7.5	6	SF08040SV	SF08040MV
8	60	6.5-7.5	6	SF08060SV	SF08060MV
8	80	6.5-7.5	6	SF08080SV	SF08080MV
8	100	6.5-7.5	6	SF08100SV	SF08100MV
8	120	6.5-7.5	6	SF08120SV	SF08120MV
8	150	6.5-7.5	6		SF08150MV
8	200	6.5-7.5	6	SF08200SV	SF08200MV
9	20	7.5-8.5	6	SF09020SV	SF09020MV
9	30	7.5-8.5	6	SF09030SV	SF09030MV
9	40	7.5-8.5	6	SF09040SV	SF09040MV
9	60	7.5-8.5	6	SF09060SV	SF09060MV
9	80	7.5-8.5	6	SF09080SV	SF09080MV
9	100	7.5-8.5	6	SF09100SV	SF09100MV
10	20	8.5-9.5	6	SF10020SV	SF10020MV
10	30	8.5-9.5	6	SF10030SV	SF10030MV
10	40	8.5-9.5	6	SF10040SV	SF10040MV
10	60	8.5-9.5	6	SF10060SV	SF10060MV
10	80	8.5-9.5	6	SF10080SV	SF10080MV
10	100	8.5-9.5	6		SF10100MV



For more information please visit www.cardinalhealth.co.uk



Important Information:

For Healthcare Professionals Only. Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information. CORDIS, Cordis LOGO and S.M.A.R.T. are trademarks of Cardinal Health and may be registered in the US and/or in other countries. © 2020 Cardinal Health. All Rights Reserved. 100551514-2 03/2020