

SLEEK™

RX PTA Balloon Dilatation Catheter

The 0.014"
rapid-exchange balloon
dedicated to infrapopliteal
lesions



Cordis™
A Cardinal Health company

SLEEK™ Catheter - rapid-exchange balloons for trackability across infrapopliteal lesions




Ultra low profile balloons designed specifically for accessing distal lesions

- 40-220 mm balloon lengths
- 2-4 mm balloon diameters
- 150 cm shaft lengths


Deliverability and crossability in infrapopliteal lesions

- Trackability in a single-operator, rapid-exchange platform
- Proprietary SLX™ lubricious coating on the balloon and the distal shaft
- Low profile 4F sheath introducer compatibility for all sizes
- Stainless steel hypotube shaft
- Nylon balloon material combines flexibility and durability with lower compliance
- High rated burst pressures of up to 16 atm
- Short balloon shoulders for precise, controlled dilatation



SLEEK™ Catheter - balloon performance for intrapopliteal lesions

SLEEK™ Catheter - superior balloon performance for intrapopliteal lesions



* Data on file. Based on bench testing on SLEEK™ Catheter and other competitors in US (4x120mm, 0.014" balloons) (internal reference: 11576423)

Product specifications

- Balloon material: Nylon
- Balloon and distal shaft coating: SLX™ (hydrophobic silicone)
- Radiopaque balloon markers: Swaged platinum iridium (2)
- Proximal shaft diameter: 2.1F
- Proximal shaft material: Stainless steel hypotube shaft with lubricious coating
- Distal shaft diameter: 3.6F (2.0-3.0 mm), 3.9F (3.5-4.0 mm)
- Distal shaft material: Pebax / Nylon
- Hub: Polycarbonate with standard luer fitting
- Maximum guidewire diameter: 0.014" (0.355 mm)

SLEEK™ Catheter - the 0.014" rapid-exchange balloon fully dedicated to infrapopliteal lesions

Ordering Information

Balloon OD x Length (mm x cm)	Rated Burst Pressure (atm)	Recommended Sheath Fit (F)	Catalogue Number	Shaft Length (cm) 150
2 X 4	16	4	425-2004	X
2 x 6	16	4	425-2006	X
2 x 12	15	4	425-2012	X
2 x 15	15	4	425-2015	X
2 x 22	15	4	425-2022	X
2.5 x 4	16	4	425-2504	X
2.5 x 6	16	4	425-2506	X
2.5 x 12	15	4	425-2512	X
2.5 x 15	15	4	425-2515	X
2.5 x 22	15	4	425-2522	X
3 x 4	16	4	425-3004	X
3 x 6	16	4	425-3006	X
3 x 12	15	4	425-3012	X
3 x 15	15	4	425-3015	X
3 x 22	15	4	425-3022	X
3.5 x 4	16	4	425-3504	X
3.5 x 6	16	4	425-3506	X
3.5 x 12	15	4	425-3512	X
3.5 x 15	15	4	425-3515	X
3.5 x 22	15	4	425-3522	X
4 x 4	16	4	425-4004	X
4 x 6	16	4	425-4006	X
4 x 12	15	4	425-4012	X
4 x 15	15	4	425-4015	X
4 x 22	15	4	425-4022	X

For more information please visit www.cardinalhealth.co.uk

